

Secure your business with the **AED7000** défibrillator

Very light (2kg), of the size of a leaf A4, in its resistant and ergonomic suitcase of transport, Starker is delivered with : 2 pairs of electrodes with simple connections, 1 adult + 1 pediatric (8-year-old children or 25 kg), self-adhesive and pre-made gel, of low impedance.

Starker allows to treat all the patients. Lithium battery supplied, guaranteed 5 years from its putting in function, deliver until 120 shocks in 150 joules.



The smallest Autonomous Defibrillator, **AED7000 v2**

Even lighter (1,5kg), smaller than a leaf A4, in its resistant and ergonomic suitcase of transport, Starker is delivered with: 2 pairs of electrodes with simple connections, 1 adult + 1 pediatric (8-year-old children or 25 kg), self-adhesive and pre-made gel, of low impedance.

Starker allows to treat all the patients. Lithium battery supplied, guaranteed 5 years from its putting in function, frees until 120 shocks in 150 joules.



Best Quality/Price report

The Defibrillator Instructor SPO2, **DM7000 Pro**

Conceived for clinical circles and hospital workers, Starker DM7000 pro, is a compact and light tool including the following main functions.

- Multi-Functional : basic configuration is very complete
- Defibrillator AED or Manuel : intuitive and easy modes of operations of uses
- Patient monitor : multi-paramétrique - SPO2 : very functional modules
- Color screen TFT of 16,5 cms,
- Printing data fast and clear on thermal paper.



STARKER DM7000 ANSWERS YOUR CLINICAL NEEDS, OF BASIC LIFE SUPPORT (BLS) TO ADVANCED CARDIAC LIFE SUPPORT (ACLS)